

# **VioTaq™** DNA Polymerase

(VT1001)

## **Description:**

Recombinant VioTaq DNA Polymerase is expressed and purified from *E. coli* to a very high homogeneity without nuclease contamination. It consists of a single polypeptide with a molecular weight of about 92 kD. VioTaq is thermostable at 95 °C with a half-time of 30 minutes and replicates DNA at 70-74 °C. This enzyme has 5'-3' exonuclease and extra A addition activity without 3'-5' exonuclease activity. It has a remarkable DNA polymerization activity with a very low error frequency (~10<sup>-5</sup>). VioTaq DNA polymerase is most suitable for large-scale PCR screening at a low cost.

## **Product Contents:**

VioTaq DNA Polymerase (stored at -20 °C), 500 U (2.5 U/μl); 10X PCR Buffer, containing 20 mM MgCl<sub>2</sub>.

## **Storage Buffer:**

20 mM Tris-HCl, pH 8.0; 0.1 mM EDTA; 1 mM DTT; 1.0% Triton X-100; 50% glycerol

## **Unit Definition:**

One unit incorporates 10 nmole of dNTP into acid-precipitable material in 30 minutes at 74 °C.

## **Basic PCR protocol:**

1. Suggested set-up of a PCR reaction

<b>Components</b>	<b>Volume</b>	<b>Final Conc.</b>
10X PCR buffer plus MgCl <sub>2</sub>	5 μl	1X (with 2.0 mM MgCl <sub>2</sub> )
10 mM each dNTP	1 μl	0.2 mM each
10 μM 5' Primer	2.5 μl	0.5 μM
10 μM 3' Primer	2.5 μl	0.5 μM
Template DNA	1-10 μl	< 0.5 μg
VioTaq (2.5 U/μl)	1 μl	2.5 U
Add sterile ddH <sub>2</sub> O to	50 μl (Final volume)	

2. Mix contents and overlay with 50 μl of mineral oil. When using a thermal cycler with a heated lid, do not need to add mineral oil.

3. Suggested PCR amplification cycles

Initial Denaturation: 94 °C for 3-5 min.

Run 25-35 cycles of the following steps:

\*Denature at 94 °C for 30-60 sec.

\*Anneal at Temp (start as 5 °C below T<sub>m</sub>) for 1 min.

\*Elongate at 72 °C for 1 min. (for 1-2 kb product)  
(Adjust elongation time based on product size, 1-2 kb per 1 min.)

Final Elongation: 72 °C for 10 min.

Maintain at 4 °C or store at -20 °C.