

User Bulletin

96 wells *miTotal*TM RNA Extraction Miniprep

System

Enhanced microRNAs Purification

Viogene[®] 96 wells *miTotal*TM RNA Extraction Miniprep System provides a high throughput method to purify total RNA from various samples such as cells from culture, tissues, whole blood, plasma, serum, biological fluids containing RNA virus, *etc.*. A simple 96 wells spin-column based method can isolate large RNAs, siRNAs, **microRNAs**, and viral RNAs without the time-consuming procedure of phenol/chloroform extraction and ethanol precipitation.

Sample Preparation Time: ~ 30 minutes, Operation under biosafety hood.

Downstream Application

- * Northern blotting
- * Poly A⁺ RNA selection
- * cDNA synthesis
- * RT-PCR
- * Transcription profiling

Product Contents

| Cat. No | VTR96 1001 |
|-------------------------------|-------------|
| Preps | 96 x 4 |
| VRX Buffer | 130ml |
| WS Buffer (RNA) | 2 x 45ml/bt |
| RNase-free ddH ₂ O | 25ml |
| VTR96 column plate | 4 |
| VTR96 Collection plate | 8 |
| Protocol | 1 |

Shipping & Storage

Viogene 96 wells *miTotal*TM RNA Extraction Miniprep System is stable at 20-25°C for one year. The kit should be stored in a dry place and kept away from direct sunlight.

Must-read Notes:

Please read the following notes before starting the procedures.

- Add 180 ml/bt of 98-100% ethanol into each WS Buffer bottle (two bottles per kit) when first open.
- All plastic wares and containers should be treated properly to make sure RNase-free. Gloves should be worn when handling RNA.
- Buffers provided in this system contain irritants. Appropriate safety apparels such as gloves and lab coat should be worn.
- **All centrifugation steps should be done at full speed of a swing-bucket centrifuge (3,000-5000 x g or 4,000-6000 rpm).**
- All procedure should be done at room temperature (20-25°C).

Protocol:

Buffer Preparation

Add 180 ml/bt of 98-100% ethanol into each WS Buffer bottle when first open.

Sample Preparation

Various guidelines are provided for monolayer cells, suspension cells, and biological fluids (whole blood, plasma, serum, buffy coat, CSF, semen, saliva, and body fluids).

Approximate confluent cell number per culture area.

| Culture Container | Well/Flask Surface | Cell Number |
|----------------------|-------------------------------|---------------------------|
| 6-well plate | 9-10 cm ² | 0.5-1x10 ⁶ |
| 12-well plate | 4 cm ² | 4-5x10 ⁵ |
| 24-well plate | 2 cm ² | 1-3x10 ⁵ |
| 96-well plate | 0.3-0.6 cm² | 4-5x10⁴ |
| T25 Culture Flask | 25 cm ² | 2-3x10 ⁶ |
| T75 Culture Flask | 75 cm ² | 0.5-1x10 ⁷ |
| T175 Culture Flask | 175 cm ² | 2-3x10 ⁷ |

Monolayer Cells

It is recommended to process 1x10⁴ – 5x10⁶ animal cells (per prep).

1. To lyse cells, add 100 µl of VRX buffer per cm² dish/plate (see table above) directly on the surface area of the adherent monolayer cell culture, mix well by pipetting, transfer the mixture into a RNase-free microcentrifuge tube or 96 plate (if use 96-well plate), and completely lyse the cells by vortexing (for tube sample) or pipetting (for 96-well plate sample).

Example: We suggest adding 100 µl VRX per well of a 96-well plate (table above).

2. Centrifuge the mixture at 3,000-5000 x g for 10 minutes to remove particulates, and then transfer the supernatant into an RNase-free microcentrifuge tube or collection plate. Proceed with RNA Purification.

Suspension Cells

It is recommended to process 1×10^4 – 5×10^6 animal cells (per prep).

1. Spin down cells by centrifugation (up to 5×10^6 animal cells). Carefully remove the supernatant, and resuspend cells in 100 μ l (for 96 wells plate sample) of VRX buffer to lyse the cell by vortexing (for tube sample) or pipetting up and down (for 96 wells plate sample).
2. Centrifuge the mixture at 3,000-5000 $\times g$ for 10 minutes to remove particulates, and then transfer the supernatant into an RNase-free microcentrifuge tube or collection plate. Proceed with RNA Purification.

Biological Fluids

100 μ l of biological liquid per prep (whole blood, plasma, serum, buffy coat, CSF, semen, saliva, biological fluids containing RNA virus and body fluids) can be processed.

1. For 100 μ l (one volume) of biological liquid, add 300 μ l (three volumes) of VRX buffer to mix the biological liquid by vortexing.
2. Centrifuge the mixture at 3,000-5000 $\times g$ for 10 minutes to remove particulates, and then transfer the supernatant into an RNase-free microcentrifuge tube or collection plate. Proceed with RNA Purification.

Note: When transfer supernatant of **whole blood** or **plasma**, avoiding bottom red blood phase

RNA Purification

1. Add **one volume isopropanol** (98-100% not provided) directly to one volume sample homogenate (1:1) in VRX buffer. Mix well by vortexing or pipetting.
2. Set 96-well miTotal RNA column plate on a collection plate.
3. Transfer the mixture into the wells of **96-well miTotal™ RNA column plate** and centrifuge at 3,000-5000 $\times g$ for 5 minutes. Discard the flow-through in the collection plate.
For biological fluids, repeat loading 400 μ l of the mixture into the wells twice.
4. **Wash** the 96 well column/collection plate **twice** with **400 μ l of WS buffer** (ethanol added) by centrifugation at 3,000-5000 $\times g$ for 5 minutes. Discard the flow-through. To complete removal of the residuals washing buffer, centrifuge the column for an **additional 10 minutes**, as residual ethanol may cause the low A_{260}/A_{230} and inhibit reverse transcriptase activity. Transfer the column carefully into an RNase-free elution plate (not provided).
Note: At this point, RNA samples can be in-column DNase treated optionally, by adding 50 μ l of DNase I cocktail (5 U DNase I, 5 μ l 10X DNase I reaction buffer, and 45 μ l RNase free water), incubate at 25-37°C for 10 minutes, after the incubation, add 100 μ l VRX/ 100 μ l isopropanol (total 200 μ l) mix into the column, centrifuge for 1 minute, discard the flow-through, then process to Step 3 and 4 of the RNA Purification.
5. Add 50 μ l of RNase-Free water **directly to the column matrix** and centrifuge at max speed for 5 minutes.
The eluted RNA can be used immediately or stored at -70 °C.

Troubleshooting

| Observation | Possible cause | Comments/suggestions |
|----------------------------|------------------------------|--|
| Low or no RNA yield | Inefficient lysis of sample | <p>Ensure that complete homogenization and disruption of samples.</p> <p>Use sufficient starting materials as suggested. Decrease it, if more than suggested amount of the sample used.</p> |
| | In sufficient sample storage | <p>Keep samples frozen until RNA extraction, avoiding samples freeze and thaw.</p> <p>If problem persist, use fresh samples and process immediately.</p> <p>Storage samples at -20°C or colder - 70°C/liquid nitrogen.</p> |
| Column clogging | | <p>Too much starting materials used. Decrease the samples to suggested sample amount.</p> |
| | | <p>Incomplete sample and VRX sample mixing. Try longer vortexing time of the mixture.</p> |
| | | <p>Check lysate for any tissues or particle remaining. Remove particles by centrifuge for 5 minutes/ 10,000xg and transfer supernatant to a new tube.</p> |
| | | <p>Ensure that complete homogenization and disruption of samples.</p> <p>Use sufficient starting materials as suggested. Decrease it, if more than suggested amount of the sample used</p> |

| | | |
|------------------------------|---------------------|---|
| | | Centrifugation temperature too low. The centrifugation temperature should be 20–25°C. |
| Degraded RNA, smear | RNase contamination | <p>Samples stored or handled incorrectly.</p> <p>Keep samples frozen until RNA extraction. Whenever possible, fresh samples should be used and processed immediately.</p> <p>Storage samples at -20°C or colder - 70°C/liquid nitrogen or RNA Stabilization reagent immediately after cells harvesting.</p> <p>Freeze & thaw cycles should be avoided.</p> <p>Wear RNase-free gloves during all procedures.</p> <p>Use only sterilized and RNase-free glasses and plasticwares.</p> |
| | | |
| DNA contamination | | Apply DNase I treatment of the suggested step. |
| No enzymatic reaction | Residues of ethanol | Before RNase-free water elution step, ensure the additional 3 minutes spin step to eliminate the residues ethanol from the washing step. |

Serial Products

| Product Name | Cat. No. | Preps |
|--|----------|-------|
| * miTotal TM RNA Miniprep | VTR1001 | 50 |
| | VTR1002 | 250 |
| * VioTotal TM Plant RNA Miniprep | PVTR1001 | 50 |
| | PVTR1002 | 250 |

* Different size of Viogene **miTotal**TM and **VioTotal**TM RNA Extraction Systems.